



Your Way To A Better Life

Introducing the finest products for Physical Therapy and Home Use!



TENS 260



ETG 255



Nerve Stimulation With Stimulation Current.

Doctors' Functions:

- Impulse width
- Impulse Frequency
- Used Programs
- Therapy Day
- Duration of Therapy
- Intensity of Therapy
- 13 Preprogramed applications
 - 5 Adjustable program applications
 - 2 Channels with adjustable intensity
- LED display for each channel
- Blue display background illumination
- Application timer function
- Program controlled frequency range 1 to 100 Hz

EMS/TENS Device ETG 255.

- Electric muscle stimulation through EMS application
- Pain alleviation through stimulation current/ TENS
- 2 channels with adjustable intensity
- LED display for each channel
- 15 user programs for muscle stimulation
- 3 user programs for nerve stimulation
- Beginner, Advanced, Professional mode
- Timer function and display background illumination



Get the devices that the professionals use.

Rest easy knowing that you purchased the #1 recommended home devices by Doctors, Physical Therapist, and Chiropractors.

CALL TODAY!

The best quality at the best price!



Greater Concepts By Design 866-867-2763
Open Monday-Friday 8am-6pm Sat 9am-1pm Sun. Closed